

DaVinci Feed Tube Set

DaVinci Feed Tube Set - 70 Series - 1 ml Catlog No.: 70-04-23-29




Figure:1

Each 1 ml DaVinci feed tube set syringe kit includes Metal Body, Luer lock nozzle, Plunger, Glass barrel, Integrated valve system, Feed Tube, Sinker valve, Feed cannula, Spare kit (Springs, Valve balls & Washers), Detailed operating manual.

Specifications:

1. Plunger tightnes without 'O' Ring.
2. Barrel made of Borosilicate glass
3. Clear visibility of Dose
4. Barrel ID 6.8 mm \pm 0.10 mm
Length 47.25 mm + 0.10 mm
5. Plunger made of Stainless Steel.
6. HS Code: 90183100




Figure:2

Table 1

Parts in Contact with Liquid For FTS			
S.No.	Component	Part Number	Material
1	Feed Tube	23	Silicone
2	Glass Barrel	7	Borosilicate glass
3	Washer, Barrel	15	FPM
4	Plunger	6	Stainless Steel
5	Valve System	11	Ni-Cr Brass
6	Washer, Valve	13	PTFE
7	Spring	12	Stainless steel
8	Ball Valve	17	Stainless steel
9	Sinker	24	Ni-Cr Brass
10	Feed cannula	26	Ni-Cr Brass

TOP SYRINGE
MANUFACTURING COMPANY PRIVATE LIMITED

10/11, Prospect Chambers Annexe, 317-21, Dr. D.N. Road, Fort, Mumbai - 400 001, Maharashtra, India.

Tel.: +91-22-6520 8849/50, +91-22-6520 6102/03 Fax: +91-22-2204 5352 • www.glass-syringe.com • sales@glass-syringe.com