

# DaVinci Feed Tube Set - 70 Series - 5 ml

## Catlog No.: 70-08-23-29




Figure:1

Each 5 ml DaVinci feed tube set syringe kit includes Metal Body, Luer lock nozzle, Plunger, Glass barrel, Integrated valve system, Feed Tube, Sinker valve, Feed cannula, Spare kit (Springs, Valve balls & Washers), Detailed operating manual.

### Specifications:

1. Plunger tightnes without 'O' Ring.
2. Barrel made of Borosilicate glass
3. Clear visibility of Dose
4. Barrel ID 15.5 mm  $\pm$  0.10 mm  
Length 47.40 mm + 0.10 mm
5. Plunger made of Stainless Steel.
6. HS Code: 90183100

Table 1

S.No.	Component	Part Number	Material
1	Feed Tube	23	Silicone
2	Glass Barrel	7	Borosilicate glass
3	Washer, Barrel	15	FPM
4	Plunger	6	Stainless Steel
5	Valve System	11	Ni-Cr Brass
6	Washer, Valve	13	PTFE
7	Spring	12	Stainless steel
8	Ball Valve	17	Stainless steel
9	Sinker	24	Ni-Cr Brass
10	Feed cannula	26	Ni-Cr Brass

#### SPARE PARTS

- |                  |                   |
|------------------|-------------------|
| 1. Piston rod    | 13. Washer        |
| 2. Sprint        | 14. Nozzle        |
| 3. Socket        | 15. Barrel washer |
| 4. Nut           | 17. Valve ball    |
| 5. Cap           | 23. Tubing        |
| 6. Piston        | 24. Sinker        |
| 7. Barrel        | (except 20 mL)    |
| 8. Casing        | 25. Vent cannula  |
| 9. Casing        | 26. Feed needle   |
| 10. Valve body   | 41. Handle        |
| 11. Intake valve | 44. Screw         |
| 12. Spring       | 46. PTFE sleeve   |
|                  | 50. Intake valve  |

Figure:2